



RÜSCH®

Care at Home



RÜSCH® MMG™ INTERMITTENT CATHETER CLOSED SYSTEM

The RüsCh® MMG™ Intermittent Catheter Closed System is clinically-proven to reduce urinary tract infections.¹ Simply remove the protective cap, load the pre-lubricated catheter into the Guardian Tip™ Introducer, insert introducer into the urethra and pass catheter into the bladder.

Available in a full range of sizes and styles, the RüsCh® MMG™ Intermittent Catheter Closed System can meet the needs of most individuals using a closed system to manage their bladder.



AT A GLANCE

- Available in kits and singles
- Full range of sizes from 6-16 Fr.
- Available in straight, coude and female lengths
- Ready to use pre-lubricated catheter
- Color-coded funnel to easily identify size
- For non-touch catheterization

Teleflex®



KIT CONTENTS

- Pre-lubricated catheter in a self-contained 1500 mL urine collection bag, except 6 Fr., 8 Fr., and 10 Fr. in a 900 mL bag
- Underpad
- One pair of gloves (not made with natural rubber latex)
- Gauze
- Three PVP swabsticks

ORDERING INFORMATION

RÜSCH® MMG™ INTERMITTENT CATHETER CLOSED SYSTEM KIT – 100/BOX

PRODUCT#	SIZE	DESCRIPTION
RLA-62-3	6 Fr.	Straight
RLA-82-3	8 Fr.	Straight
RLA-102-3	10 Fr.	Straight
RLA-122-3	12 Fr.	Straight
RLA-122-3C	12 Fr.	Coudé
RLA-122-3O	12 Fr.	Soft
RLA-142-3	14 Fr.	Straight
RLA-142-3C	14 Fr.	Coudé
RLA-142-3F	14 Fr.	Female Length
RLA-142-3O	14 Fr.	Soft
RLA-162-3	16 Fr.	Straight
RLA-162-3C	16 Fr.	Coudé
RLA-182-3	18 Fr.	Straight
SONK-141-3*	14 Fr.	Straight
SONK-142-3**	14 Fr.	Straight
The following kits contain Red Rubber Catheters (latex)		
RLA-122-3R	12 Fr.	Straight
RLA-142-3R	14 Fr.	Straight
RLA-142-3RC	14 Fr.	Coudé
RLA-162-3R	16 Fr.	Straight

*This kit contains only gauze and PVP swabsticks

** This kit contains only gauze, PVP swabsticks and gloves (not made with natural rubber latex)

TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990

REFERENCE

1. Bennett CJ, Young, MN et al. The Effect of Urethral Introducer Tip Catheters On The Incidence of Urinary Tract Infection Outcomes in Spinal Cord Injured Patients, *Journal of Urology*, Vol. 168 No. 2 August 1997.

Rx only

Teleflex, Guardian Tip, MMG and Rüsch are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.
© 2015 Teleflex Incorporated. All rights reserved. MC-000068

RÜSCH® MMG™ INTERMITTENT CATHETER CLOSED SYSTEM SINGLES – 100/BOX

PRODUCT#	SIZE	DESCRIPTION
ONC-8	8 Fr.	Straight
ONC-10	10 Fr.	Straight
ONC-12	12 Fr.	Straight
ONC-12C	12 Fr.	Coudé
ONC-14	14 Fr.	Straight-Soft PVC
ONC-14C	14 Fr.	Coudé
ONC-16	16 Fr.	Straight
ONC-16C	16 Fr.	Coudé
SONC-14	14 Fr.	Straight-Standard PVC
SONC-14F	14 Fr.	Female Length

The following products contain Red Rubber Catheters (latex)

RONC-14	14 Fr.	Straight Red Rubber
RONC-16	16 Fr.	Straight Red Rubber

Contains pre-lubricated catheter in self contained collection bag without supplies.